



RSID: Semen & Saliva

+ Reduced Incubation Time

~ 10 Seconds using Vortex mixer



To determine if RSID-Semen and RSID-Saliva are compatible with shorter sample extraction times, a series of time course experiments were undertaken with control swabs, aged samples (several years old), trace semen samples, and semen on fabrics/ trace saliva samples, and salivas on fabrics.

These data clearly demonstrated that similar results could be obtained from all tested sample types using incubation times as short as 10 seconds (*with shaking*) to as long as 1 hour

(to view the data, go to <http://galantos.eu/rsid-semen.html> and <http://galantos.eu/rsid-saliva.html>)

Room temperature extraction of forensic samples for a minimum of 10 seconds *with shaking* is sufficient for detection of semenogelin with RSID-Semen and α -amylase RSID-Saliva. Longer incubation times (i.e., 5-60 minutes) are optional.



+ Required Protocol for RSID-Semen



Maximum sensitivity of RSID™-Semen can be obtained with the following protocol:

Extraction volume: 300-400 μ L of RSID buffer

Incubation time: minimum 10 seconds with shaking, longer times are optional

Use up to 20 μ L of the extract plus 80 μ L of RSID-Semen buffer to a final of 100 μ L.

Apply to cassette.

Read / Score results at 10 minutes.

+ Stability and Storage

RSID cassettes can be shipped and stored at room temperature.

The RSID products can be shipped and stored at room temperature.

Once the bottles with the buffers have been opened they should be stored in the fridge.

If possible store the buffers in the fridge.

© 2018 Galantos Genetics GmbH