

# SPERM HY-LITER™ SOS

## Solution Staining Protocol



### Kit Provided Solutions:

SAMPLE PREPARATION Solution  
BLOCKING Solution  
SPERM HEAD STAINING Solution  
MOUNTING Media

yellow bottle cap (*addition of DTT required before use*)  
red bottle cap  
green bottle cap  
blue bottle cap

---

### User-prepared solutions:

**0.1 M Dithiothreitol (DTT) Solution** – Use Independent Forensics 1 M DTT or Prepare 1 M DTT according to the provided DTT solution recipe below. Dilute 1 M DTT stock 1:10 in DiH<sub>2</sub>O (*i.e., 10 µL 1M DTT + 90 µL DiH<sub>2</sub>O*).

#### 1 M DTT Solution Recipe

Final Volume:	10 mL	100 mL
DTT	1.54 g	15.4 g
H <sub>2</sub> O	~7 mL	~70 mL
1 M KOH	1.1 mL	11 mL

Confirm that the prepared DTT solution is ~pH 8.0. Add DiH<sub>2</sub>O to the final volume, aliquot and freeze. Aliquots may be frozen and thawed twice before discarding.

#### Sample Preparation Solution (SPS)+ DTT - Prepare SPS + DTT daily before use:

Prepare sufficient SPS + DTT master mix for all of your samples, plus 1 extra. Each sample will require 40 µL of sample prep solution; dispense sufficient SP into a tube and remove [(#samples + 1) x 40] µL into a fresh tube.

Add 5 µL of 0.1M stock DTT solution per (sample +1) to complete the SPS + DTT master mix.

For example, if three (3) samples are to be stained, add 4-5 drops Sample Preparation Solution to a tube and dispense 160 µL of SP into a fresh tube. Add 20 µL of 0.1 M DTT to the measured SP solution.

*Sperm cells in solution are significantly more sensitive to DTT relative to sperm stained on slides. Accurate concentrations and incubation times are particularly important in the SOS technique.*

---

### Staining Procedure:

Use the following protocol only with **SPERM HY-LITER™ SOS** solutions; Do not substitute reagents from the Express kit

#### Preparation of extracts:

- Extract swab or fabric cutting in 1.7 mL microcentrifuge tube in 0.5 mL of PBS. Incubate on bench top for 1-2 hrs with occasional mixing.
- Recover cells using Spin basket: transfer the soaked swab head or fabric cutting into a spin basket, place the basket with the substrate into its same 1.7 mL tube. Centrifuge for 5 minutes at max rpm (~13,000 x g).

Carefully remove supernatant from cell pellet – to avoid disturbing the pellet, leave behind about 50  $\mu\text{L}$ , do not remove all the supernatant. **For staining, use 10-20  $\mu\text{L}$  of resuspended cell pellet for staining in solution.**

*Suggested control swabs:*

Negative – buccal swab pellet resuspended in 50-100  $\mu\text{L}$  of PBS

Positive – semen cell pellet resuspended in 50  $\mu\text{L}$  of PBS.

- Prepare three controls samples for staining:

A. Buccal only (10  $\mu\text{L}$  of resuspended Buccal cell pellet in a tube)

B. Sperm only (10  $\mu\text{L}$  of resuspended sperm cell pellet in a tube)

C. Buccal + sperm (10  $\mu\text{L}$  each of resuspended buccal and sperm cell pellets).

**A. Prepare 0.1 M DTT and Sample Preparation Solution + DTT: DTT stock for protocol on slides is 1 M. The DTT stock for in-solution staining is 0.1 M and must be diluted from the 1 M stock.**

### **B. Steps for Staining In Solution:**

1. **Sample Preparation:** Add 45  $\mu\text{L}$  of SPS/DTT solution (prepared as above) to each sample to be stained. Incubate at room temperature for 5 min with occasional gentle mixing.
2. **Block:** Add 1 drop **BLOCKING Solution (red bottle cap)** to each sample. Incubate at room temperature for 5 minutes with occasional mixing.
3. Centrifuge tube for 3 minutes at 7,000 RPM and gently remove most of supernatant (all but ~50  $\mu\text{L}$ ) being careful not to disturb the cell pellet.
4. **Stain:** Add 1 drop **SPERM HEAD STAINING Solution (green bottle cap)** to each sample. Resuspend cell pellet by pipetting up and down and incubate at room temperature for 5 minutes.
5. Centrifuge for 3 minutes at 7,000 RPM and gently remove most of supernatant (all but ~50  $\mu\text{L}$ ), being careful not to disturb the cell pellet.
6. Resuspend pellet in 10-20  $\mu\text{L}$  PBS and deposit entire volume onto 11 mm sample window of a microscope slide. Allow to completely dry before proceeding to next step.
7. **Mount:** Add one drop of **MOUNTING Media (blue bottle cap)** to each sample window. Gently place provided cover slip over each sample window. Place slide between a folded paper towel and gently press down to position cover slip and remove excess mounting media. Mounting media will semi-harden after 20 minutes at room temperature. Slides may be stabilized by outlining the cover slip with clear nail polish. Slides may be visualized immediately and are stable for days.

### **C. Visualize:**

Stained slides must be visualized using a fluorescence microscope fitted with appropriate filters. Cell nuclei, including epithelial and sperm, can be visualized using DAPI-compatible filters. Human sperm heads can be visualized using fluorescein or Alexa 488 compatible filters. Slides may be scanned using 10x, 20x, 40x or 100x objectives at the operator's discretion.